# AlaxoStent


DA	BRUGSANVISNING	01
DE	GEBRAUCHSANWEISUNG	04
EN	INSTRUCTIONS FOR USE	07
ES	MANUAL DE INSTRUCCIONES	10
FI	KÄYTTÖOHJEET	13
FR	MANUEL D'UTILISATION	16
GR	ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ	19
HR	UPUTE ZA UPOTREBU	22
HU	HASZNÁLATI ÚTMUTATÓ	25
IT	ISTRUZIONI PER L'USO	28
NL	GEBRUIKSAANWIJZING	31
NO	BRUKERVEILEDNING	34
PT	INSTRUÇÕES DE USO	37
SV	BRUKSANVISNING	40
TR	KULLANMA TALIMATLARI	43

**REF**

AP004 AlaxoStent C

[service@alaxo.com](mailto:service@alaxo.com)

[www.alaxo.com](http://www.alaxo.com)


1

- DA Komponenter af AlaxoStent  
 DE Komponenten des AlaxoStent  
 EN Components of AlaxoStent  
 ES Componentes del AlaxoStent  
 FI AlaxoStent-tuotteen komponentit  
 FR Composants de l'AlaxoStent  
 GR Συστατικά μέρη της AlaxoStent  
 HR Komponente proizvoda AlaxoStent  
 HU Az AlaxoStent komponensei  
 IT Componenti dello AlaxoStent  
 NL Onderdelen van de AlaxoStent  
 NO Komponentene til AlaxoStent  
 PT Componentes de AlaxoStent  
 SV AlaxoStent-systemets delar  
 TR AlaxoStent Bileşenleri


2


- DA Sådan sættes de to dele af klipselukningen sammen  
 DE ineinanderschieben des Clipverschlusses  
 EN Fitting into each other of the clip lock  
 ES Acople del cierre de clip  
 FI Pihitlukon osien sovittaminen toisiinsa  
 FR Télescopage des deux parties du clip de fermeture  
 GR Σύνδεση των δύο μερών με το κλιπ  
 HR Umetanje dijelova kopče za spajanje jedan u drugi  
 HU a klipzár összekapcsolása  
 IT Collegare la chiusura a clip  
 NL In elkaar schuiven van de clipsluitering  
 NO Delene av klemlåsen settes i hverandre  
 PT Encaixando-se mutuamente no clipe de bloqueio  
 SV Passa ihop klicklåsets ändrar  
 TR Klipsli kilidin birbirine takılması

3

- DA Nitinolgitteret indsættes i indføringsrøret  
 DE Einführen des Nitinolgeflechts in den Einführschlauch  
 EN Insertion of the nitinol braid into the introduction tube  
 ES Inserción de la malla de nitinol en el cateter de inserción  
 FI Nitinolipunoksen vieminen sisäänvientiputkeen  
 FR Introduction du treillis de nitinol dans le tube d'introduction  
 GR Εισαγωγή του πλέγματος νιτινόλης στο σωληνάκι εισαγωγής  
 HR Uvođenje nitinolske mrežice u cijev vodilicu  
 HU a nitinol szövet behúzása a bevezető tömlőbe  
 IT Inserire la rete in nitinol nel tubo di inserzione  
 NL Invoeren van het nitinolvlechtwerk in de invoerslang  
 NO Stikk nitinolnettet inn i innføringsslangen  
 PT Inserção do entrançado nitinol no interior do tubo introdução  
 SV För in nitinolnätet helt i införingsslangen  
 TR Nitinol ağ örgünün giriş tüpüne sokulması

4

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## EN INSTRUCTIONS FOR USE

### 1 Instructions

The AlaxoStent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

#### 1.1 CE Mark

The product bears the CE mark. The conformity evaluation assures that the product complies with the requirements of the EU directive for medical devices 93/42/EWG. According to the classification rules in Annex IX of the EU directive the AlaxoStent is a Class I product.

#### 1.2 General Instructions

The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product.

In case of damage resulting from use of third-party components or consumables no warranty claims apply.

Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.

The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.

Full or partial reproduction subject to prior written approval by Alaxo.

#### 1.3 General safety instructions

The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.

It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.

The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.

The product must not be accessible to children.

The user has to inspect the product for functional safety and good condition each time prior to application.

The user has to be familiar with handling of the product.


The cleaning instructions have to be observed and obeyed. Other cleaning procedures must not be applied. The cleaning spray (Prontosan® Wound Spray) must only be used for disinfectant cleaning of the AlaxoStent and the accessories. Do not swallow!

### 1.4 Symbols and indications

One Prontosan® Wound Spray bottle may be used for 12 months after first use.

The AlaxoStent may be used for max. 18 months after the first use.

In the instructions for use, on packagings and on the product the following figures and symbols are used:

 Consult instructions for use	 Batch code	 Order number
 Temp. limitation: store between 4° and 50°C	 Manufacturing date: Month/Year	 Use by: Month/Year
 CE Mark	 Manufacturer	 Store inaccessible for children

### 1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoStent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoStent in the throat and its removal.

The attending doctor has to advise on potential risks and complications connected to the AlaxoStent therapy and has to evaluate if the AlaxoStent therapy is applicable to the specific patient.

### 2 Product information

#### 2.1 Intended use

The AlaxoStent is a medical device and self-applied by the patient before going to sleep. The AlaxoStent is introduced via a nasal passage into the throat and positioned there in order to keep open the airway during sleep by mechanical splinting.

The AlaxoStent is a prescription medical device. The attending doctor has to be sufficiently qualified in treatment of sleep apnea.

Only one AlaxoStent is introduced through one nasal passage, never two or more at the same time.

The AlaxoStent is clinically tested and CE marked for treatment of patients with obstructive sleep apnea or snoring.

The attending doctor has to determine the disease state of the patient and if the AlaxoStent may be a useful therapy for this patient. The physician further has to examine if the nasal passages of the patient allow introduction of the AlaxoStent. The physician shall inform the patient which of the two nasal passages (left or right) is preferred for introduction of the AlaxoStent.

In case of known lung diseases a pneumologic medical examination is required to evaluate if the AlaxoStent therapy can be applied to the specific patient.

The AlaxoStent must be used only by one single patient.

The AlaxoStent is suitable for repeated use up to 18 hours per application.

The AlaxoStent must be used only by patients of age of 18 years and higher.

The AlaxoStent is not suitable for life-supporting application in emergency care.

## 2.2 Contraindications

- Inappropriate structure of the nasal passages, especially if leading to problems during insertion of the introduction tube into a nasal passage
- Strong deviations of the nasal septum
- Dry rhinitis
- Disposition for nose hemorrhage
- Empty Nose Syndrome
- Shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes

## 2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Persons with allergic reactions to polyhexanide or undecylenamidopropyl-betaine may suffer an allergic response to the cleaning spray.
- Do not use the AlaxoStent if this leads to pain! Do not terminate therapy without consulting your attending doctor.

## 2.4 Unintended use

- The only permitted use is splinting of the airway by introduction through a nasal passage into the throat. Other uses are prohibited.
- The AlaxoStent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoStent braid must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

## 2.5 Product description

The AlaxoStent comprises the following components (see  1):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3]
- 4) Cleaning tube [4]
- 5) Fixation tape
- 6) 0.1% Polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray)


## 3 Scope of Delivery

- |   |                        |
|---|------------------------|
| 1 Nitinol braid                                       | 1 Pusher               |
| 1 Introduction tube                                   | 1 Cleaning tube        |
| 1 Fixation tape                                       | 1 Instructions for use |
| 1 Bottle cleaning spray (Prontosan® Wound Spray 75ml) |                        |



## 4 Use

The AlaxoStent is a medical device and has been developed and clinically tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed in order to maintain the nitinol braid in good order over the allowed period of use of max. 18 months.

## 4.1 Preparation of the AlaxoStent for introduction

1. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see  2).
2. The other end of the pusher [2b] is inserted into black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube and the end of the stent shortly before the black end. (see  3)



## 4.2 Introduction of the AlaxoStent into the throat

1. The introduction tube is carefully introduced with the black end ahead into the right or left nostril (tilt the head back!) and advanced horizontally through the nasal passage according to  4 down to the throat. Exact positioning of the nitinol braid is explained to the patient by the prescribing doctor. In each case the AlaxoStent stays away from the epiglottis.
2. The smooth material and the bent form of the introduction tube serves for good steerability at the end of the nasal passage down to the throat and is optimally adapted to the human anatomy. Nevertheless, differences in simplicity of the introduction process may arise e.g. from anatomic variability from patient to patient.
3. In case of dry mucosal tissue in the nose pushing forward of the introduction tube may lead to increased friction. To wetten the nasal mucosa, after consultation of the attending physician, a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be applied horizontally into the nasal passage. The introduction tube smoothly moves forward along the wet mucosal tissue with low friction.
4. After reaching the right position the pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands first in the throat and then in the nasal passage and attaches itself towards the walls of the throat.
5. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid. The AlaxoStent thus is in its position of use (see  4). Cleaning of the introduction tube is described in Section 4.4.
6. To secure the position of the AlaxoStent in the throat the part of the nitinol braid protruding from the nose is affixed with the fixation tape at the upper lip, the nose or the cheek (see  5). It has to be made sure that the fixation tape securely adheres to the nitinol braid.
7. Application of the AlaxoStent is completed and the patient can go to sleep.
8. Due to the opening force of the nitinol braid the throat is kept open and a collapse of the airway is prevented. Differences in the efficacy of splinting of the airway may be observed from patient to patient.

## 4.3 Removal of the AlaxoStent

After sleep the fixation tape is cautiously removed and the nitinol braid retracted from the nose and throat.

## 4.4 Cleaning of the AlaxoStent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm (not hot !!!) flowing drinking water until all potentially adhering mucus is removed.
2. The outer surface of the introduction tube is sprayed with the cleaning spray (Prontosan® Wound Spray) which is also sprayed into the end of the introduction tube. The cleaning spray has to act for one minute and then the introduction tube is let dry.
3. The **nitinol braid** is rinsed in lukewarm (not hot !!!) flowing drinking water immediately after retraction from the nose until all potentially adhering mucus is mostly removed.
4. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed over one end of the tube (do not pull!!!). Subsequently push the whole nitinol braid over the tube as shown in  6. Removal of contaminations from the nitinol braid with the fingers has to be done only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube) as otherwise the sensitive nitinol braid may be damaged.
5. Once the nitinol braid is completely cleaned it is cautiously pushed off the cleaning tube with the thumb and one finger without exerting any pressure from the lower end to the upper end (see  7). Please do not pull the closed end nor the clip lock of the nitinol braid as it will contract due to its design and cannot be pulled off the cleaning tube!
6. Subsequently the nitinol braid is sprayed with the cleaning spray (Prontosan® Wound Spray) which has to act for one minute. Then the nitinol braid is rinsed once more in lukewarm flowing drinking water, again sprayed with the cleaning spray and let dry.
7. The **cleaning tube** and the **pusher** are rinsed, too, in lukewarm flowing drinking water, sprayed with the cleaning spray (Prontosan® Wound Spray) which has to act for one minute and then let dry.

The three-step cleaning procedure for the AlaxoStent serves for disinfection. By extensive rinsing the contaminations are efficiently removed. The 0.1% polyhexanide / 0.1% betaine cleaning spray (Prontosan® Wound Spray) additionally reaches all interspaces due to its modified surface tension. The ingredients possess excellent mucosa compatibility and a broad pathogen spectrum.

## 4.5 Storage of the AlaxoStent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The AlaxoStent shall not be stored in direct sunlight and at temperatures below 4°C or above 50°C.

## 4.6 Period of use

1. The AlaxoStent after first use may be used for max. 18 months for nightly application.
2. After this period an increased risk of breakage of a nitinol wire due to material wear exists. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. Inspection is easiest performed on the cleaning tube during the cleaning procedure.
4. A cleaning spray bottle (Prontosan® Wound Spray) may be used for 12 months after first use. Expiry date see print on label.
5. After expiry of these allowed periods of use a new AlaxoStent has to be employed.

## 5 Disposal

The AlaxoStent can be disposed of in domestic waste.

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Alaxo GmbH