# Alaxo Hybrid Stent

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EN: Instructions for Use

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
FR: Mode d'emploi

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AP011 Alaxo Hybrid Stent

[www.alaxo.com](http://www.alaxo.com)

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## 1. Instructions

The Alaxo Hybrid Stent is a medical device.

### 1.1 CE Mark

The product bears the CE mark. The declaration of conformity assures that the product complies with the requirements of the Regulation (EU) 2017/745 (MDR). According to the classification rules in Annex VIII of the MDR it is a Class I product.

### 1.2 General Instructions

The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product. In case of damage resulting from use of third-party components or consumables no warranty claims apply.

Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.

The instructions for use corresponds to the design of the product and the applicable safety standards as of going to print. All rights are reserved for the specified parts and components, processes, names and products.

Full or partial reproduction subject to prior written approval by Alaxo.

### 1.3 General safety instructions

The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.

It is not allowed to modify or reconstruct the product. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.

The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary Alaxo does not take any liability for damages resulting from insufficient packaging.

The product must not be accessible to children.

The user has to inspect the product for functional safety and good condition each time prior to application. The user has to be familiar with handling of the product.

The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied.

## 1.4 Symbols and information

The Alaxo Hybrid Stent set may be used for max. 18 months after first use.

In the instructions for use, on the packaging and on the product the following figures and symbols are used:

Order number

Manufacturer

Lot number

Manufacturing date

Expiry date

Consult Instructions for Use

Store inaccessible for children

Temperature limitation:  
store between lower and upper temperature

Medical device

Unique device identifier

## 1.5 Precautionary measures

Correct handling of the Alaxo Hybrid Stent is described in this instructions for use and has to be adhered to.

Potential risks and complications connected to the use of the Alaxo Hybrid Stent are described in this instructions for use.

## 2 Product information

### 2.1 Intended use

- The Alaxo Hybrid Stent set is a medical device and self-applied by the user before going to sleep. The long stent (Alaxo Hybrid Stent) is introduced through one nasal passage into the throat. The short stent (AlaxoLito Plus Nasal Stent) is inserted into the other nasal passage, but not into the throat.
- The Alaxo Hybrid Stent set is intended for splinting of the nasal passage and the upper pharynx.

- The AlaxoHybrid Stent set must be used only by a single person.
- The Alaxo Hybrid Stent set is suitable for repeated use up to 18 hours per application.
- The Alaxo Hybrid Stent set must be used only by persons of age of 18 years and higher.
- The Alaxo Hybrid Stent set is not suitable for life-supporting application in emergency care.

## 2.2 Contraindications

- Strong deviations of the nasal septum
- Disposition for nose hemorrhage
- Empty Nose Syndrome
- Shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes

## 2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Do not use the product if this leads to pain. See a doctor for medical examination.

## 2.4 Unintended use

- The only permitted use is splinting of the airway by introduction into the nasal passage and the throat. Other uses are prohibited.
- The stents must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The stents must not be pulled over any body parts or limbs.
- The stents may only be used for the approved medical purposes. They must not be used for any other purposes, particularly non-medical purposes.

## 2.5 Product description

The Alaxo Hybrid Stent set comprises the following components (see Fig. 1):

A: Alaxo Hybrid Stent, B: AlaxoLito Plus Stent

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices) = stent
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3] with black end [3a]
- 4) Cleaning tube [4]

### 3 Scope of delivery

- 1 Alaxo Hybrid Stent nitinol stent
- 1 Alaxo Hybrid Stent introduction tube
- 1 Alaxo Hybrid Stent pusher
- 1 Alaxo Hybrid Stent cleaning tube
- 1 AlaxoLito Plus nitinol stent
- 1 AlaxoLito Plus introduction tube
- 1 AlaxoLito Plus pusher
- 1 AlaxoLito Plus cleaning tube
- 1 Instructions for use

### 4 Use

The Alaxo Hybrid Stent set is a medical device for self-application by the user. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed to in order to maintain the nitinol braid in good order during the allowed period of use of max. 18 months.

#### 4.1 Preparation of the stents for introduction

1. The stent is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see Fig. 2: Fitting into each other of the clip lock).
2. The other end of the pusher [2b] is inserted into the black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the stent is carefully and completely pulled into the introduction tube using the pusher whereby the braid self-compresses. The clip lock and the end of the stent shall be located in the introduction tube just before the black end (see Fig. 3: Insertion of the nitinol braid into the introduction tube).

#### 4.2 Insertion of the stents

1. Insertion of the Alaxo Hybrid Stent into the throat: The introduction tube is carefully introduced with the black end ahead into one nostril and advanced horizontally in direction of the cheek bone through the nasal passage according to Fig. 4 (Scheme for positioning of the Alaxo Hybrid Stent; A=Nasal passage, B=Soft palate, C=Tongue base) as far as the clip lock still being outside of the nose (tilt the head back to facilitate insertion).
2. Insertion of the AlaxoLito Plus Nasal Stent into the nose: The introduction tube is carefully introduced with the black end ahead into the other nostril and advanced horizontally in direction of the cheek bone into the nasal passage as far as the clip lock still being outside of the nose (see Fig. 5: Scheme for positioning of the AlaxoLito Plus Nasal Stent).

3. In case of dry mucosa in the nose pushing forward of the introduction tube may lead to increased friction. To wetten the nasal mucosa a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid (but no deswelling drug) may be applied horizontally into the nasal passage. The introduction tube smoothly moves forward along the wet mucosal tissue with low friction.
4. After reaching the right position the pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands and attaches to the mucosa.
5. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the stent. The respective stent thus is in its operating position as determined by the ball-shaped section being placed at the nasal alar.
6. The Alaxo Hybrid Stent splints the upper throat and thereby prevents a collapse of the airway in the region of the soft palate. Further, it splints the nasal passage and therewith optimizes nasal breathing.
7. Nasal breathing on the other side of the nose is optimized by the AlaxoLito Plus Nasal Stent
8. Differences in the efficacy of splinting of the airway may be observed for different users.

#### **4.3 Removal of the stents**

For removal the two stents are cautiously pulled out of the nose. Should the nasal mucosa have dried out it is recommended to first moisten the mucosa. For that purpose a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be sprayed into the nose.

#### **4.4 Cleaning of the stents**

1. After retraction from the nose the introduction tube is rinsed in lukewarm (not hot!) flowing drinking water until all potentially adhering mucus is removed.
2. The stent is rinsed in lukewarm (not hot!) flowing drinking water immediately after retraction from the nose until all potentially adhering mucus is mostly removed. In case of dried mucus on the stent immersion in lukewarm drinking water for 15 min may be used to dissolve the mucus and then continue with rinsing in flowing drinking water.
3. The stent is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed under lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed over one end of the tube (do not pull!). Subsequently push the nitinol braid over the tube as far as shown as shown in Fig. 6 (Pushing of the nitinol braid over the cleaning tube; left: Alaxo Hybrid Stent; right:

AlaxoLito Plus Nasal Stent). Removal of contaminations from the stent with the fingers has to be done only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube) and cautiously to prevent damage to the sensitive stent braid.

4. Once the stent is completely cleaned it is cautiously pushed off the cleaning tube with the thumb and one finger without exerting any pressure from the lower end to the upper end (see Fig. 7: Pushing off of the nitinol braid from the cleaning tube). Do not pull the closed end nor the clip lock of the nitinol braid as it will contract due to its design and cannot be pulled off the cleaning tube!
5. The cleaning tube and the pusher shall be rinsed, too, in lukewarm flowing drinking water.

#### **4.5 Storage**

1. The stent shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The Alaxo Hybrid Stent set shall be stored not in direct sunlight and not at a temperature below 4°C and above 50°C.

#### **4.6 Period of use**

1. The Alaxo Hybrid Stent set after first use may be used for max. 18 months for nightly application.
2. After this period an increased risk of breakage of a nitinol wire due to material fatigue exists. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. Inspection is easiest performed on the cleaning tube during the cleaning procedure.

#### **5 Disposal**

The Alaxo Hybrid Stent set can be disposed of in domestic waste.

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Alaxo GmbH