
AlaxoLito

BRUGSANVISNING
GEBRAUCHSANWEISUNG
INSTRUCTIONS FOR USE
MANUAL DE INSTRUCCIONES
KÄYTTÖOHJEET
MANUEL D'UTILISATION
ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ
UPUTE ZA UPOTREBU
HASZNÁLATI ÚTMUTATÓ
ISTRUZIONI PER L'USO
GEBRUIKSAANWIJZING
BRUKERVEILEDNING
INSTRUÇÕES DE USO
BRUKSANVISNING
KULLANMA TALIMATLARI

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

REF

AP001 AlaxoLito
AP002 AlaxoLito Plus
AP005 AlaxoLito Xtreme
AP006 AlaxoLito Sport

service@alaxo.com
www.alaxo.com

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EN INSTRUCTIONS FOR USE

1 Instructions

The AlaxoLito Nasal Stent is a CE marked prescription medical device. It must only be used according to the briefing of the attending doctor.

1.1 CE Mark

The product bears the CE mark. The declaration of conformity assures that the product complies with the requirements of the Regulation (EU) 2017/745 (MDR). According to the classification rules in Annex VIII of the MDR the AlaxoLito Nasal Stent is a Class I product.

1.2 General Instructions

- The instructions for use is a component of the product. It has to be accessible to the user near to the product. Close adherence to these instructions is precondition for intended use and correct handling of the product.
- In case of damage resulting from use of third-party components or consumables no warranty claims apply.
- Alaxo is only liable for the product in respect of safety, reliability and functionality if it is used in compliance with the instructions for use.
- The instructions for use corresponds to the design of the product and the applicable safety norms as of going to print. All rights are reserved for the specified parts and components, processes, names and products.
- Full or partial reproduction subject to prior written approval by Alaxo.

1.3 General safety instructions

- The product has been developed and designed by Alaxo in such a manner that hazards during intended use are eliminated as far as possible. Nevertheless, we strongly recommend to strictly adhere to the following safety precautions in order to further reduce any residual risks to a minimum.
- It is not allowed to modify or reconstruct the product. In that case the CE mark lapses. Use of altered products is prohibited. In addition, Alaxo cannot take any warranty or liability for modified or reconstructed products. In order to assure safe application of the product users are liable for adherence to the instructions.
- The original packaging should be stored for any potential return shipment. Only the original packaging ensures optimal protection of the product during transport. Should return shipment become necessary during the warranty period Alaxo does not take any liability for damages resulting from insufficient packaging.
- The product must not be accessible to children.
- The user has to inspect the product for functional safety and good condition each time prior to application.
- The user has to be familiar with handling of the product.
- The cleaning instructions have to be observed and obeyed to. Other cleaning procedures must not be applied.

DA

DE

EN

ES

FI

FR

GR

HR

HU

IT

NL

NO

PT

SV

TR

1.4 Symbols and indications

- The AlaxoLito Nasal Stent may be used for 18 months after first use.
- In the instructions for use, on packagings and on the product the figures and symbols in  are used.

1.5 Precautionary measures

The attending doctor has to explain handling of the AlaxoLito Nasal Stent to the patient. Especially, the doctor has to explain correct application and positioning of the AlaxoLito stent in the nose and its removal.

The attending doctor has to advise on potential risks and complications connected to the AlaxoLito Nasal Stent therapy and has to evaluate if the AlaxoLito therapy is applicable to the specific patient.

2 Product information

2.1 Intended use

- The AlaxoLito Nasal Stent is a medical device and self-applied by the patient. The AlaxoLito stent is positioned in the nasal cavity in order to mechanically splint the nasal airway.
- The AlaxoLito Nasal Stent splints the anterior middle or inferior nasal passage. The AlaxoLito Plus/Xtreme/Sport Nasal Stent splints the middle or inferior nasal passage over its entire length. The AlaxoLito Nasal Stent can be used to treat nasal valve or nasal alar collapse and turbinate hyperplasia as well as potentially resulting snoring.
- The AlaxoLito Nasal Stent is a CE marked and prescription medical device.
- The attending doctor has to examine if the AlaxoLito Nasal Stent is an appropriate therapy for the patient and if the nasal cavities and passages of the patient are appropriate for application of the AlaxoLito stent.
- The AlaxoLito Nasal Stent must be used only by one single patient.
- The AlaxoLito Nasal Stent is suitable for repeated use up to 18 hours per application.
- The AlaxoLito Nasal Stent must be used only by patients of age of 18 years and higher.
- The AlaxoLito Nasal Stent is not suitable for life-supporting application in emergency care.

2.2 Contraindications

- Inappropriate structure of the nasal passages, especially if leading to problems during insertion of the introduction tube into a nasal passage

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

- Strong deviations of the nasal septum
- Dry rhinitis
- Disposition for nose hemorrhage
- Empty Nose Syndrome
- Shall not be used by persons with reduced mental ability, with reduced motor skills, with psychic restrictions, mental disorders or spastic syndromes

2.3 Warnings

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this stent.
- Do not use the AlaxoLito Nasal Stent if this leads to pain! Do not terminate therapy without consulting your attending doctor.

2.4 Unintended use

- The only permitted use is splinting of the airway in the nasal cavity. Other uses are prohibited.
- The AlaxoLito Nasal Stent must not be introduced into orifices of the body other than the nose (e.g. ears or mouth).
- The AlaxoLito stent must not be pulled over any body parts or limbs.
- The product may only be used for the approved medical purposes. It must not be used for any other purposes, especially non-medical purposes.

2.5 Product description

The AlaxoLito Nasal Stent comprises the following components (see  1):

- 1) Self-expanding nitinol braid [1] with clip lock [1a] (nitinol is a shape-memory metal, a nickel-titanium alloy widely used for medical devices)
- 2) Pusher [2] with clip lock (the end with the two holes) [2a]
- 3) Introduction tube [3] with black end [3a]
- 4) Cleaning tube [4]

3 Scope of Delivery

- 2 Nitinol braids
- 1 Introduction tube
- 1 Pusher
- 1 Cleaning tube
- 1 Instructions for use

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

4 Use

The AlaxoLito Nasal Stent is a medical device and has been developed and tested for self-application by the patient. The nitinol braid consists of thin and sensitive metal wires so that the braid needs to be treated carefully and gently. The following instructions for handling and cleaning have to be obeyed to in order to maintain the nitinol braids in good order over the allowed period of use of 18 months.

DA

DE

EN

ES

FI

FR

GR

HR

HU

IT

NL

NO

PT

SV

TR

4.1 Preparation of the AlaxoLito stent for introduction

1. The nitinol braid is connected to the pusher [2] by fitting into each other the two parts [1a] and [2a] of the clip lock (see  2).
2. The other end of the pusher is inserted into the black end [3a] of the introduction tube [3] and the pusher is moved through the introduction tube. Then the nitinol braid is carefully and completely pulled into the introduction tube using the pusher whereby it self-compresses. The clip lock shall be located in the introduction tube. (see  3)

4.2 Introduction of the AlaxoLito stent into the nose

1. The introduction tube containing the nitinol braid is carefully introduced with the black end ahead into the left or right nostril as far as the clip lock still being outside of the nose ( 4).
2. The pusher is held with the one hand and the introduction tube is slowly retracted over the pusher with the other hand. The nitinol braid self-expands and thereby attaches itself towards the walls of the nose. Thus the AlaxoLito stent is safely and stably positioned.
3. Once the introduction tube is fully removed the clip lock is pulled apart and thereby the pusher is separated from the nitinol braid.
4. Subsequently the other stent is introduced into the other nostril in the same way.
5. Cleaning of the introduction tube is described in Section 4.4.

4.3 Removal of the AlaxoLito stent

1. For removal the two nitinol braids are cautiously pulled out of the nose. Should the nasal mucosa have dried out it is recommended to first moisten the mucosa.
2. For that purpose a nasal spray containing sea salt, dexpanthenol and/or hyaluronic acid may be sprayed into the nose.

4.4 Cleaning of the AlaxoLito Nasal Stent

1. After retraction from the nose the **introduction tube** is rinsed in lukewarm

- (not hot !!!) flowing drinking water until all potentially adhering mucus is removed.
2. The **nitinol braid** is rinsed in lukewarm (not hot !!!) flowing drinking water immediately after pulling out of the nose until all potentially adhering mucus is mostly removed.
 3. The nitinol braid is pushed over the cleaning tube [4] and potentially adhering contaminations are carefully removed in lukewarm flowing drinking water with the fingers. For bringing the nitinol braid onto the cleaning tube its open end is cautiously pushed - obliquely fitted - over one end of the tube (do not pull!!!). Subsequently push the whole nitinol braid over the tube as shown in  5. Removal of contaminations from the nitinol braid with the fingers has to be done only in longitudinal direction, never in cross direction (i.e. not around the cleaning tube) as otherwise the sensitive nitinol braid may be damaged.
 4. Once the nitinol braid is completely cleaned it is cautiously pulled off the cleaning tube.
 5. The **cleaning tube** is rinsed, too, in lukewarm flowing drinking water and let dry.
 6. Should the **pusher** be contaminated it is also cleaned by rinsing in lukewarm flowing drinking water, spraying with the cleaning spray which has to act for one minute and let dry.

4.5 Storage of the AlaxoLito Nasal Stent

1. The dried nitinol braid shall not be stored in the introduction tube to prevent it from unnecessary pressure load.
2. The AlaxoLito shall not be stored in direct sunlight and at temperatures below 4°C or above 50°C.

DA
DE
EN
ES
FI
FR
GR
HR
HU
IT
NL
NO
PT
SV
TR

4.6 Period of use

1. The AlaxoLito after first use may be used for 18 months for one application per day.
2. After this period an increased risk of breakage of a nitinol wire due to material fatigue exists, especially for the round wires at the open end of the nitinol braid. Therefore, usage is no longer allowed.
3. The nitinol braid must be inspected for any damage prior to each use. If breakage of a wire is observed the nitinol braid must not be used anymore. Inspection is easiest performed on the cleaning tube during the cleaning procedure.

5 Disposal

The AlaxoLito can be disposed of in domestic waste.

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Alaxo GmbH